

Keysight Technologies

5600LS AFM High-resolution Imaging Molecular-level Understanding of n-Alkanes Self-Assembly onto Graphite

Application Note




Figure 1. AFM topographic image of n-C₃₆H₇₄ on graphite. Scan size: 350 nm × 350 nm.




Figure 2. AFM topographic image of n-C₃₆H₇₄ on graphite. Scan size: 55 nm × 55 nm.

The adsorption of organic molecules from solution onto a solid surface has attracted tremendous attention as it is fundamentally associated to many phenomena of both industrial and academic relevance. Governed by an intricate balance between adsorbate-substrate and adsorbate-adsorbate interactions, spontaneous self-organization of molecules at the interface may lead to the formation of delicate ultrathin films with nanometerscale ordered surface structures due to the molecular-level packing in a particular way. Various techniques have been reported to investigate those organic layers. For instance, differential scanning calorimetry (DSC) provides an effective means to probe the surface phase behavior. The structural information normal to the interface, to some extent, can be extracted from neutron reflectivity measurements. So far, the real time 3D structural characterization with atomic or molecular scale details of the assemblies especially the top layer comes mainly from scanning probe microscopy techniques. In this application brief, the

capability of atomic force microscopy (AFM) to directly visualize soft thin film materials with true molecular resolution is demonstrated using self-assembly of n-C₃₆H₇₄ molecules on graphite as an example. All the data displayed here are acquired from a Keysight Technologies, Inc. 5600LS system with a 90 μm large scanner. A typical AFM topographic image of n-C₃₆H₇₄ upon adsorption on HOPG is shown in Figure 1, from which rich information about this sample is revealed. First, molecules are aligned on the substrate with long-range order and exhibit a striped morphology. Second, existence of local defect areas is captured. As can be seen, the whole image is divided into four segments by two long and one short domain boundaries and a protrusion island is observed near the bottom location. Those surface features correspond to adsorbed molecules in an amorphous state. Third, two different orientations of the molecular packing are identified. The stripes in the left-upper corner are exactly 60° rotated with respect to those in the remaining three domains, reflecting

the impact of underneath substrate (i.e., the 6-fold symmetry graphite) on n-C₃₆H₇₄ alignment on the surface. Figure 2 is another n-C₃₆H₇₄ / graphite topography image with a larger magnification to deliver important information at single molecular level. It shows that the stripe width is about 4.5 nm, which is matching well with the molecular length of n-C₃₆H₇₄ with a fully extended configuration. Furthermore, the linear backbones (i.e., hydrocarbon chains with an all-tans configuration) of individual n-C₃₆H₇₄ molecules are resolved at lower part of the image. These results unambiguously illustrate that n-C₃₆H₇₄ molecules are lying down and orderly aligned on graphite to form a lamellar packing structure.

In conclusion, AFM is a powerful surface characterization tool with an unprecedented high-resolution. Materials surface structures with sub-5 nm size in lateral dimensions can be resolved clearly, thus making it possible to achieve molecular-level understanding of the adsorption behavior of long-chain molecules at the solid-solution interface.

Jing-jiang Yu, Ph.D.
Keysight Technologies

AFM Instrumentation from Keysight Technologies

Keysight offers high-precision, modular AFM solutions for research, industry, and education. Exceptional worldwide support is provided by experienced application scientists and technical service personnel. Keysight's leading-edge R&D laboratories are dedicated to the timely introduction and optimization of innovative and easy-to-use AFM technologies.

For more information on Keysight Technologies' products, applications or services, please contact your local Keysight office. The complete list is available at: www.keysight.com/find/contactus

Americas

Canada	(877) 894 4414
Brazil	55 11 3351 7010
Mexico	001 800 254 2440
United States	(800) 829 4444

Asia Pacific

Australia	1 800 629 485
China	800 810 0189
Hong Kong	800 938 693
India	1 800 112 929
Japan	0120 (421) 345
Korea	080 769 0800
Malaysia	1 800 888 848
Singapore	1 800 375 8100
Taiwan	0800 047 866
Other AP Countries	(65) 6375 8100

Europe & Middle East

Austria	0800 001122
Belgium	0800 58580
Finland	0800 523252
France	0805 980333
Germany	0800 6270999
Ireland	1800 832700
Israel	1 809 343051
Italy	800 599100
Luxembourg	+32 800 58580
Netherlands	0800 0233200
Russia	8800 5009286
Spain	800 000154
Sweden	0200 882255
Switzerland	0800 805353
	Opt. 1 (DE)
	Opt. 2 (FR)
	Opt. 3 (IT)
United Kingdom	0800 0260637

For other unlisted countries:
www.keysight.com/find/contactus
(BP-09-23-14)

www.keysight.com/find/afm



This information is subject to change without notice.
© Keysight Technologies, 2009 - 2014
Published in USA, July 31, 2014
5990-3956EN
www.keysight.com